



*Your Satisfaction is our Technology*



**High Frequency**

**VECTOR Series...3.5/5/6/6R**

## TECHNICAL SPECIFICATION

Specification	High Frequency Mobile/Fixed X-Ray Generator			
Model	VECTOR-3.5	VECTOR-5	VECTOR-6	VECTOR-6/R
Max Output(Kw)	3.5 Kw	5Kw	6Kw	6Kw(Rotating)
Kv Range	40-110Kv	40-110Kv	40-110Kv	40-110Kv
mA Range	20-80mA	20-100mA	20-120mA	20-150mA
Focal Spot	0.6/1.5		0.3/0.6	
mAs Range	0.2 to 250mAs			
Time	4ms to 5 Sec.			
X-Ray Tube	Stationary Anode dual focus			Rotating Anode dual focus
Digital Display	Kv,mA&mAs			
Control Panel	Membrane Key Board,LCD			
APR Facility	Available			
Collimator	Manually/Motorized Light Beam Diaphragm(LBD)-Optional			
Power Requirement	1 Phase 200/220/230 VAC with 15Amp Power Socket			
Tube Stands	Counter balanced Mobile/Spring balanced Mobile/Ceiling to wall Stand			
Battery/Capacitor Bank	Optional			

The company reserved the rights to add /delete/Modify in details of specification and appearance mentioned here.

Quality and Safety Certifications.....



**ALLANS**  
MEDICAL SYSTEMS  
PVT. LTD.

No.-1,Thakur Lane,Opp. Datta Mandir  
Behind Disha Hospital, Aslpha Village  
Ghatkopar(W), Mumbai-400084  
Tel:-022-25128274 / 8879780695  
Email:- [allansmedisys@gmail.com](mailto:allansmedisys@gmail.com)  
Web:-[www.allansmedisys.in](http://www.allansmedisys.in)

Represented by

